

pN-Z-TK₂ 2-step, nlacZ targeting vector

Polylinker

3'-CGAACCATGGGCCCCTA
ATCGATAAGCTAGCTTGGTACCCGGGGAT
 ClaI NheI* KpnI* BamHI

GGAGAT-5'
CCTCTAGAGTCGACGGATCCGGGGAATTC
 XbaI Sall* BamHI EcoRI

CCCAGTCTCAGGATCCACCATGG
 BamHI NcoI

* Note: AvrII, XbaI, and SpeI have ends that are cohesive with NheI
 Asp718 and Acc65 also recognize this sequence;
 BsiW1 & BsrG1 are cohesive with Acc65
 XhoI is cohesive with Sall

The ATG of the NcoI site is the start of nuclear LacZ

An anti-sense PCR primer (~70° annealing temperature) is shown in italics #582

Suggestions: Clone promoter fragment of convenient PCR length (~1 kb) into polylinker in correct orientation. Clone largest piece possible into Xho, Spe, Not sites. Plan to linearize construct before electroporation with SstII or AscI (if your DNA has both of these sites, convert them to something unique first). If you don't want the lacZ gene, clone one piece of homologous DNA between NheI, KpnI, Sall and HindIII sites.

All restriction sites shown on the map are unique

