

FIB-SEM / SEM Array Tomography Sample Preparation of Cell Monolayer

Note: This protocol is designed to enhance signal for backscatter electron (BSE) imaging of epoxy embedded cell monolayer at low accelerating voltages (1-5 keV).

Day 1					
Steps		Time	Temp	Check	
Fixation (I)	2.5% GA in 0.1M CB	5 m x 1	36.5 °C	<input type="checkbox"/>	
	EMS, Cat #15960	25 m x 1	RT	<input type="checkbox"/>	
Washing	0.1M CB	3 m x 3	RT	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
Fixation (II)	1% K ₄ Fe(CN) ₆ + 1% OsO ₄ in 0.1M CB	15 m x 1	ICE	<input type="checkbox"/>	
	0.2ml 10% potassium ferrocyanide K ₄ Fe(CN) ₆ , 1%				
	0.5ml 4% OsO ₄ , 1%				
	0.5ml 0.4M CB, 0.1M				
	0.8ml MQW				
	2 ml Total				
Washing	Mili-Q/Distilled Water	3 m x 3	RT	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
TCH	1% Thiocarbohydrazide (TCH) in Mili-Q/Distilled Water	5 m x 1	RT	<input type="checkbox"/>	
	<ul style="list-style-type: none"> ✓ Make fresh and filter right before use! ✓ Dissolve 0.02g TCH with 2ml Mili-Q/Distilled Water in glass vial at 60 °C! ✓ A toxic compound made by the reaction of carbon disulfide with hydrazine (hydrazinolysis). ✓ It is used in the silver proteinate specific staining of carbohydrates. 				
	Washing	Mili-Q/Distilled Water	3 m x 3	RT	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	OsO ₄	1% OsO ₄ in Mili-Q/Distilled Water	15 m x 1	RT	<input type="checkbox"/>
Washing	Mili-Q/Distilled Water	3 m x 3	RT	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
en bloc Stain	1% UA in Mili-Q/Distilled Water	Overnight	4 °C	<input type="checkbox"/>	
Day 2					
Washing	Mili-Q/Distilled Water	3 m x 3	RT	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
Walton's Lead Asparate Staining		10 m x 1	60 °C	<input type="checkbox"/>	
	<ul style="list-style-type: none"> ✓ Prepare an aspartic acid stock solution by dissolving 0.998 g of L-aspartic acid in 250 ml of MQW. ✓ Dissolve 0.05g L-aspartic acid in 12.5ml boiled distilled water. ✓ Note: the aspartic acid will dissolve more quickly if the pH raised to 3.8. ✓ This stock solution is stable for 1-2 months if refrigerated. ✓ To make a lead asparate solution dissolve 0.066 g of lead nitrate in 10ml of aspartic acid stock and pH adjusted to 5.5 with 1N KOH/NaOH. ✓ 12 drops of 1N KOH/NaOH using 1 ml measurer disposable pipette. ✓ Note: the lead asparate solution is placed in a 60 °C oven for 30 minutes. 				
Washing	MQW	3 m x 3	RT	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	

FIB-SEM / SEM Array Tomography Sample Preparation of Cell Monolayer

Dehydration	10% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	30% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	50% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	70% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	80% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	90% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	95% EtOH	5 m x 1	ICE	<input type="checkbox"/>
	100% EtOH (I)	5 m x 1	RT	<input type="checkbox"/>
	100% EtOH (II)	5 m x 1	RT	<input type="checkbox"/>
Day 3				
Polymerization	Absolute Durcupan	30 m x 1	RT	<input type="checkbox"/>
	Durcupan : 100% EtOH = 1 : 1	1 h x 1	RT	<input type="checkbox"/>
	Durcupan : 100% EtOH = 2 : 1	4 h x 1	RT	<input type="checkbox"/>
	Absolute Durcupan	Overnight	RT	<input type="checkbox"/>
Day 4				
Re-sizing	<ul style="list-style-type: none"> ✓ ~ 1.2 mm³ of resin blocks ✓ Attach on resin cylinder for SEM Array Tomography ✓ Attach on SEM pin stub for FIB-SEM imaging 			<input type="checkbox"/>
Sectioning	<ul style="list-style-type: none"> ✓ 750 x 750 x 70 nm for TEM imaging ✓ 800 x 800 x 50 nm for SEM Array Tomography 			<input type="checkbox"/>
Carbon Evaporation	<ul style="list-style-type: none"> ✓ Sections on a Si wafer 	Carbon Thread	7 nm	<input type="checkbox"/>
Sputtering	<ul style="list-style-type: none"> ✓ ~ 1.2 mm³ of resin block on SEM pin stub 	Au	10 nm	<input type="checkbox"/>

FIB-SEM / SEM Array Tomography Sample Preparation of Cell Monolayer